



TYPICAL DRAWER

NOTES

1. STANDARD ELECTRICAL SERVICE:
115/1/60, 3-WIRE CORD PLUG
208/1/60, 230/1/60,
OTHER ELECTRICAL CHARACTERISTICS
AVAILABLE
2. MAX. OPERATING TEMPERATURE:
350 DEG. F
3. AIR FLOW: HORIZONTAL, SIDE TO
SIDE
4. INTERIOR: 304 STAINLESS STEEL, WITH
CARBON STEEL BLOWER WHEEL
5. EXTERIOR: 304 BRUSHED STAINLESS STEEL



OPTIONAL LEG STAND

This drawing is the property of THE GRIEVE CORPORATION, and may not be used or reproduced without the express permission of THE GRIEVE CORPORATION

REV.	DESCRIPTION	BY	DATE
DESCRIPTION			
CUSTOMER			
MODEL DR-350 FORCED CONVECTION BENCH OVENS			
SHOP ORDER			
DRAWN BY PJD		DWG. NO.	
DATE		PRINT DATE 5/31/2013	
		DR-350	